




GAGE7 rabbit pAb

Catalog No	YP-Ab-11230
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GAGE7 GAGE12I GAGE7B
Immunogen	Synthesized peptide derived from human GAGE7 AA range: 41-91
Specificity	This antibody detects endogenous levels of GAGE7 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AL4;Cancer/testis antigen 4.7;CT4.7;G antigen 12F;G antigen 7;GAGE7;GAGE7 GAGE12I GAGE7B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in some prostate cancer tissues but not in normal prostate tissue.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night