



TRI66 rabbit pAb

Catalog No	YP-Ab-11235
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRIM66 C11orf29 KIAA0298
Immunogen	Synthesized peptide derived from human TRI66 AA range: 608-658
Specificity	This antibody detects endogenous levels of TRI66 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRI66;TRIM66 C11orf29 KIAA0298
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus. Forms discrete foci within the centromeric chromocenter and surrounding nucleoplasm. .
Tissue Specificity	
Function	function:May function as transcription repressor; The repressive effects are mediated, at least in part, by recruitment of deacetylase activity. May play a role as negative regulator of postmeiotic genes acting through CBX3 complex formation and centromere association.,similarity:Contains 1 bromo domain.,similarity:Contains 1 PHD-type zinc finger.,similarity:Contains 2 B box-type zinc fingers.,subcellular location:Forms discrete foci within the centromeric chromocenter and surrounding nucleoplasm.,subunit:Can form homo- and heterodimers. Interacts with CBX5, CBX1 and CBX3 via PxVxL motif.,
Background	


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4°over night