




ER6L2 rabbit pAb

Catalog No	YP-Ab-11240
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ERCC6L2 C9orf102
Immunogen	Synthesized peptide derived from human ER6L2 AA range: 193-243
Specificity	This antibody detects endogenous levels of ER6L2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ER6L2;ERCC6L2 C9orf102
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,mitochondrion,microtubule organizing center,
Tissue Specificity	
Function	sequence caution:Translation N-terminally extended.,similarity:Belongs to the SNF2/RAD54 helicase family.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,
Background	This gene encodes a member of the Snf2 family of helicase-like proteins. The encoded protein may play a role in DNA repair and mitochondrial function. Mutations in this gene have been associated with bone marrow failure syndrome 2. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Apr 2014],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night