




PKHA3 rabbit pAb


Catalog No	YP-Ab-11241
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	PLEKHA3 FAPP1
Immunogen	Synthesized peptide derived from human PKHA3 AA range: 189-239
Specificity	This antibody detects endogenous levels of PKHA3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHA3;PLEKHA3 FAPP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Localizes to ER-trans Golgi network (TGN) contacts sites. .
Tissue Specificity	Widely expressed.
Function	function: Binds specifically to phosphatidylinositol-4-phosphate (PtdIns4P), but not to other phosphoinositides.,similarity: Contains 1 PH domain.,tissue specificity: Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human uterus. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).