




ABCAB rabbit pAb

Catalog No	YP-Ab-11254
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ABCA11P ABCA11
Immunogen	Synthesized peptide derived from human ABCAB AA range: 45-95
Specificity	This antibody detects endogenous levels of ABCAB at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABCAB;ABCA11P ABCA11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of membrane,
Tissue Specificity	
Function	caution:In contrast to other members of the family, it is much shorter and lacks ABC transporter domains. Defined as a pseudogene by HGNC. However, proteomic data on a specific peptide from PubMed:17924679 suggest the existence of this protein.,similarity:Belongs to the ABC transporter family. ABCA subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night