




PCOTH rabbit pAb

Catalog No	YP-Ab-11263
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	PCOTH
Immunogen	Synthesized peptide derived from human PCOTH AA range: 14-64
Specificity	This antibody detects endogenous levels of PCOTH at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCOTH;prostate collagen triple helix;Protein PCOTH
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in prostate and testis. Weakly or not expressed in other tissues. Overexpressed in prostate cancers.
Function	function:May be involved in growth and survival of prostate cancer cells through the TAF-Ibeta pathway.,tissue specificity:Expressed in prostate and testis. Weakly or not expressed in other tissues. Overexpressed in prostate cancers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA, pH 9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight). 3, Secondary antibody was diluted at 1:200 (room temperature, 45 min).