




TMCC3 rabbit pAb

Catalog No	YP-Ab-11264
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMCC3 KIAA1145
Immunogen	Synthesized peptide derived from human TMCC3 AA range: 33-83
Specificity	This antibody detects endogenous levels of TMCC3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1145;Q9ULS5;Transmembrane and coiled-coil domain protein 3;TMCC3; TMCC3 KIAA1145
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Concentrates in discrete patches along peripheral endoplasmic reticulum tubules. .
Tissue Specificity	Widely expressed, with highest levels in brain, spinal cord and testis.
Function	similarity:Belongs to the TEX28 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night