




T120B rabbit pAb

Catalog No	YP-Ab-11268
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM120B
Immunogen	Synthesized peptide derived from human T120B AA range: 71-121
Specificity	This antibody detects endogenous levels of T120B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T120B;TMEM120B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus inner membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the TMEM120 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night