



ANAG rabbit pAb

Catalog No	YP-Ab-11281
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NAGLU UFHSD1
Immunogen	Synthesized peptide derived from human ANAG AA range: 146-196
Specificity	This antibody detects endogenous levels of ANAG at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANAG;NAGLU UFHSD1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Lysosome.
Tissue Specificity	Liver, ovary, peripheral blood leukocytes, testis, prostate, spleen, colon, lung, placenta and kidney.
Function	catalytic activity:Hydrolysis of terminal non-reducing N-acetyl-D-glucosamine residues in N-acetyl-alpha-D-glucosaminides.,caution:A MPS3B mutation at position 100 was erroneously reported (PubMed:9950362) as an amino acid change from Arg to His. The right amino acid change is from His to Arg., disease:Defects in NAGLU are the cause of mucopolysaccharidosis type 3B (MPS3B) [MIM:252920]; also known as Sanfilippo syndrome B. MPS3B is a form of mucopolysaccharidosis type 3, an autosomal recessive lysosomal storage disease due to impaired degradation of heparan sulfate. MPS3 is characterized by severe central nervous system degeneration, but only mild somatic disease. Onset of clinical features usually occurs between 2 and 6 years; severe neurologic degeneration occurs in most patients between 6 and 10 years of age, and death occurs typically during the second or third decade of



life.,functio

Background

This gene encodes an enzyme that degrades heparan sulfate by hydrolysis of terminal N-acetyl-D-glucosamine residues in N-acetyl-alpha-D-glucosaminides. Defects in this gene are the cause of mucopolysaccharidosis type IIIB (MPS-IIIB), also known as Sanfilippo syndrome B. This disease is characterized by the lysosomal accumulation and urinary excretion of heparan sulfate. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4° over night