




MESD rabbit pAb


Catalog No	YP-Ab-11283
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; IHC
Gene Name	MESDC2 KIAA0081 MESD UNQ1911/PRO4369
Immunogen	Synthesized peptide derived from human MESD AA range: 125-175
Specificity	This antibody detects endogenous levels of MESD at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BOCA; KIAA0081; LDLR chaperone MESD; LRP chaperone MESD; MESD; MESDC2 KIAA0081 MESD UNQ1911/PRO4369
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum . Released from apoptotic cells and shed photoreceptor outer segments. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).