




CK080 rabbit pAb

Catalog No	YP-Ab-11285
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C11orf80
Immunogen	Synthesized peptide derived from human CK080 AA range: 93-143
Specificity	This antibody detects endogenous levels of CK080 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CK080;C11orf80
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Chromosome . Localizes to meiotic chromosomes. .
Tissue Specificity	Detected in lung, spleen,colon and in skeletal muscle. Expressed in the ovaries, Fallopian tubes and uterus (PubMed:30388401).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night