




MOSC1 rabbit pAb


Catalog No	YP-Ab-11302
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;ELISA;IHC
Gene Name	MARC1 MOSC1
Immunogen	Synthesized peptide derived from human MOSC1 AA range: 156-206
Specificity	This antibody detects endogenous levels of MOSC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MARC1;MOSC1;MTARC1;MARC1 MOSC1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion outer membrane ; Single-pass type II membrane protein . Membrane ; Lipid-anchor . Mitochondrial import is mediated by AA 1-40 and requires ATP.
Tissue Specificity	
Function	cofactor:Molybdenum.,function:Probable oxidoreductase.,similarity:Contains 1 MOSC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).