




CCDC7 rabbit pAb

Catalog No	YP-Ab-11315
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CCDC7 BIOT2
Immunogen	Synthesized peptide derived from human CCDC7 AA range: 380-430
Specificity	This antibody detects endogenous levels of CCDC7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein BIOT2;Coiled-coil domain-containing protein 7;CCDC7;BIOT2;CCDC7 BIOT2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in epithelium of normal cervix and cervical cancer. Overexpressed in early and interim cervical cancer.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4°over night