




TM242 rabbit pAb

Catalog No	YP-Ab-11328
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM242 C6orf35 BM-033 HSPC287
Immunogen	Synthesized peptide derived from human TM242 AA range: 53-103
Specificity	This antibody detects endogenous levels of TM242 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM242;TMEM242 C6orf35 BM-033 HSPC287
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the UPF0463 family.,similarity:Contains 28 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night