




C1QT6 rabbit pAb

Catalog No	YP-Ab-11338
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	C1QTNF6 CTRP6 UNQ581/PRO1151
Immunogen	Synthesized peptide derived from human C1QT6 AA range: 31-81
Specificity	This antibody detects endogenous levels of C1QT6 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1QT6; C1QTNF6 CTRP6 UNQ581/PRO1151
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity: Contains 1 C1q domain., similarity: Contains 1 collagen-like domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4°over night