




PDXD2 rabbit pAb

Catalog No	YP-Ab-11345
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PDXDC2P PDXDC2
Immunogen	Synthesized peptide derived from human PDXD2 AA range: 144-194
Specificity	This antibody detects endogenous levels of PDXD2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDXD2;PDXDC2P PDXDC2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	cofactor:Pyridoxal phosphate.,similarity:Belongs to the group II decarboxylase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night