




DPPA2 rabbit pAb

Catalog No	YP-Ab-11353
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DPPA2 PESCRG1
Immunogen	Synthesized peptide derived from human DPPA2 AA range: 138-188
Specificity	This antibody detects endogenous levels of DPPA2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECSA;Pluripotent embryonic stem cell-related gene 1 protein;PESCRG1; Developmental;DPPA2;DPPA2 PESCRG1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues.
Function	function:May play a role in maintaining cell pluripotentiality.,similarity:Contains 1 SAP domain.,tissue specificity:Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night