




GPBL1 rabbit pAb

Catalog No	YP-Ab-11359
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	GPBP1L1 SP192
Immunogen	Synthesized peptide derived from human GPBL1 AA range: 170-220
Specificity	This antibody detects endogenous levels of GPBL1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Dp 1;DP1;DRTF1;DRTF1 polypeptide 1;E2F dimerization partner 1;TFDP1;transcription;GPBL1;GPBP1L1 SP192
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Possible transcription factor.,similarity:Belongs to the vasculin family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night