




FOXD2 rabbit pAb

Catalog No	YP-Ab-11364
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FOXD2 FKHL17 FREAC9
Immunogen	Synthesized peptide derived from human FOXD2 AA range: 177-227
Specificity	This antibody detects endogenous levels of FOXD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKHL17;Forkhead box protein D2;Forkhead-related protein FKHL17;Forkhead-related transcription;FOXD2;FOXD2 FKHL17 FREAC9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Kidney specific.
Function	function:Probable transcription factor involved in embryogenesis and somatogenesis.,similarity:Contains 1 fork-head DNA-binding domain.,tissue specificity:Kidney specific.,
Background	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night