




L12R1 rabbit pAb


Catalog No	YP-Ab-11373
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	LOH12CR1
Immunogen	Synthesized peptide derived from human L12R1 AA range: 9-59
Specificity	This antibody detects endogenous levels of L12R1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L12R1;LOH12CR1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Lysosome membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	Ubiquitously expressed (at protein level).
Function	similarity:Belongs to the LOH12CR1 family.,tissue specificity:Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).