




TT39C rabbit pAb

Catalog No	YP-Ab-11374
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TTC39C C18orf17
Immunogen	Synthesized peptide derived from human TT39C AA range: 229-279
Specificity	This antibody detects endogenous levels of TT39C at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HsT2697;C18orf17;TT39C;TTC39C C18orf17
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the TTC39 family.,similarity:Contains 3 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night