



TBL3 rabbit pAb

Catalog No	YP-Ab-11389
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TBL3 SAZD
Immunogen	Synthesized peptide derived from human TBL3 AA range: 486-536
Specificity	This antibody detects endogenous levels of TBL3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SAZD; TBL3; transducin (beta) like 3; Transducin beta like protein 3; TBL3 SAZD
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	similarity: Contains 13 WD repeats.,
Background	The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD40 group is a large family of proteins, which appear to have a regulatory function. It is believed that the WD40 repeats mediate protein-protein interactions and members of the family are involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation. This gene has multiple polyadenylation sites. It might have multiple alternatively spliced transcript variants but the variants have not been fully described yet. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night