




RWD2B rabbit pAb

Catalog No	YP-Ab-11395
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RWDD2B C21orf6 GL011
Immunogen	Synthesized peptide derived from human RWD2B AA range: 228-278
Specificity	This antibody detects endogenous levels of RWD2B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RWD2B;RWDD2B C21orf6 GL011
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	similarity:Contains 1 RWD domain.,tissue specificity:Ubiquitous.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night