




CAHM3 rabbit pAb

Catalog No	YP-Ab-11397
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CALHM3 FAM26A
Immunogen	Synthesized peptide derived from human CAHM3 AA range: 191-241
Specificity	This antibody detects endogenous levels of CAHM3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAHM3;CALHM3 FAM26A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Basolateral cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Specifically expressed in circumvallate taste bud cells.
Function	similarity:Belongs to the FAM26 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night