




NCK5L rabbit pAb

Catalog No	YP-Ab-11416
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NCKAP5L KIAA1602 FP1193
Immunogen	Synthesized peptide derived from human NCK5L AA range: 261-311
Specificity	This antibody detects endogenous levels of NCK5L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cep169;FP1193;KIAA1602;Nck associated protein 5 like;NCKAP5L;NCK5L; NCKAP5L KIAA1602 FP1193
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Localizes to microtubule plus ends (PubMed:26485573, PubMed:26482847). Associates with centrosomes during interphase, but dissociates from these structures from the onset of mitosis (PubMed:26549230). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night