




TM230 rabbit pAb

Catalog No	YP-Ab-11420
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM230 C20orf30 HSPC274 UNQ2432/PRO4992
Immunogen	Synthesized peptide derived from human TM230 AA range: 32-82
Specificity	This antibody detects endogenous levels of TM230 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf30; HSPC274; TMEM230; UNQ2432/PRO4992; TM230; TMEM230 C20orf30 HSPC274 UNQ2432/PRO4992
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . Early endosome . Recycling endosome . Late endosome . Cytoplasmic vesicle, autophagosome .
Tissue Specificity	
Function	similarity: Belongs to the UPF0414 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4°over night