




HECD2 rabbit pAb

Catalog No	YP-Ab-11424
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	HECTD2
Immunogen	Synthesized peptide derived from human HECD2 AA range: 505-555
Specificity	This antibody detects endogenous levels of HECD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HECT domain containing 2;HEC;HECD2;HECTD2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,cytoplasm,cytosol,
Tissue Specificity	
Function	function:Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night