




GTSF1 rabbit pAb

Catalog No	YP-Ab-11441
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	GTSF1 FAM112B
Immunogen	Synthesized peptide derived from human GTSF1 AA range: 41-91
Specificity	This antibody detects endogenous levels of GTSF1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cue11;FAM112B;gametocyte specific factor 1;GTSF1;Protein FAM112B;GTSF1 FAM112B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	similarity:Belongs to the UPF0224 (FAM112) family.,similarity:Contains 1 CHHC-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night