




U17L5 rabbit pAb

Catalog No	YP-Ab-11457
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	USP17L5
Immunogen	Synthesized peptide derived from human U17L5 AA range: 308-358
Specificity	This antibody detects endogenous levels of U17L5 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U17L5;USP17L5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Endoplasmic reticulum .
Tissue Specificity	
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., caution:The RS447 megasatellite DNA is a highly polymorphic conserved tandem repetitive sequence which contains a copy of USP17. It is present with an interindividual variation in copy number and between 20 to 103 copies can be found in the genome both on chromosome 4 and chromosome 8., function:Functions in cell apoptosis. Cleaves ubiquitin fusion protein substrates., similarity:Belongs to the peptidase C19 family. USP17 subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night