




S7A13 rabbit pAb

Catalog No	YP-Ab-11461
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	SLC7A13 AGT1 XAT2
Immunogen	Synthesized peptide derived from human S7A13 AA: 20-50
Specificity	This antibody detects endogenous levels of S7A13 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	XAT2;X amino acid transporter 2;AGT1;AGT 1;S7A13;SLC7A13 AGT1 XAT2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in the kidney.
Function	function:Mediates the transport L-aspartate and L-glutamate in a sodium-independent manner.,similarity:Belongs to the amino acid-polyamine-organocation (APC) superfamily.,tissue specificity:Expressed in the kidney.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4°over night