




FND3A rabbit pAb

Catalog No	YP-Ab-11476
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FNDC3A FNDC3 HUGO KIAA0970
Immunogen	Synthesized peptide derived from human FND3A AA range: 184-234
Specificity	This antibody detects endogenous levels of FND3A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fibronectin type-III domain-containing protein 3A;FNDC3;HUGO;FND3A; FNDC3A FNDC3 HUGO KIAA0970
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in the odontoblast and nerves in the dental pulp. Also expressed in trachea and to a lesser extent in the brain, liver, lung and kidney.
Function	similarity:Contains 9 fibronectin type-III domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night