




F100A rabbit pAb

Catalog No	YP-Ab-11480
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FAM100A PP11303
Immunogen	Synthesized peptide derived from human F100A AA range: 100-150
Specificity	This antibody detects endogenous levels of F100A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F100A; FAM100A PP11303
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity: Belongs to the FAM100 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night