




AVPI1 rabbit pAb

Catalog No	YP-Ab-11491
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	AVPI1 NPD013 PP5395
Immunogen	Synthesized peptide derived from human AVPI1 AA range: 76-126
Specificity	This antibody detects endogenous levels of AVPI1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Arginine vasopressin-induced protein 1; AVP induced protein 1; AVP-induced protein 1; NPD13; AVPI1; AVPI1 NPD013 PP5395
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function: May be involved in MAP kinase activation, epithelial sodium channel (ENaC) down-regulation and cell cycling.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4°over night