




IG2R rabbit pAb

Catalog No	YP-Ab-11505
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GIG44 PP9974
Immunogen	Synthesized peptide derived from human IG2R AA range: 21-71
Specificity	This antibody detects endogenous levels of IG2R at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IG2R;GIG44 PP9974
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in fetal and adult liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4°over night