




BWR1B rabbit pAb

Catalog No	YP-Ab-11515
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SLC22A18AS BWR1B BWSCR1B ORCTL2S SLC22A1LS
Immunogen	Synthesized peptide derived from human BWR1B AA range: 44-94
Specificity	This antibody detects endogenous levels of BWR1B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BWR1B;BWSCR1B;ORCTL2S;p27 BWR1B;SLC22A18AS;SLC22A1LS; SLC22A18AS BWR1B BWSCR1B ORCTL2S SLC22A1LS
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Most abundantly expressed in gastrointestinal tissues. Expressed at lower levels in kidney and placenta. Expressed in fetal brain, liver, placenta, kidney and lung.
Function	tissue specificity:Most abundantly expressed in gastrointestinal tissues. Expressed at lower levels in kidney and placenta. Expressed in fetal brain, liver, placenta, kidneys and lung.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4°over night