




KBTB2 rabbit pAb


Catalog No	YP-Ab-11516
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC
Gene Name	KBTBD2 BKLHD1 KIAA1489 CGI-73
Immunogen	Synthesized peptide derived from human KBTB2 AA range: 359-409
Specificity	This antibody detects endogenous levels of KBTB2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BKLHD1;KBTBD2;KIAA1489;KBTB2;KBTBD2 BKLHD1 KIAA1489 CGI-73
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cul3-RING ubiquitin ligase complex,
Tissue Specificity	
Function	similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 Kelch repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).