




TRML4 rabbit pAb

Catalog No	YP-Ab-11520
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TREML4 TLT4 UNQ9425/PRO34675
Immunogen	Synthesized peptide derived from human TRML4 AA range: 140-190
Specificity	This antibody detects endogenous levels of TRML4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TLT 4;TLT4;Trem like transcript 4 protein;TREML4;TRML4;TREML4 TLT4 UNQ9425/PRO34675
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night