




TMCC1 rabbit pAb

Catalog No	YP-Ab-11527
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMCC1 KIAA0779
Immunogen	Synthesized peptide derived from human TMCC1 AA range: 275-325
Specificity	This antibody detects endogenous levels of TMCC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0779;TMCC 1;Transmembrane and coiled-coil domains protein 1;TMCC1; TMCC1 KIAA0779
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Specifically localizes to contact sites between the endoplasmic reticulum and endosomes that are spatially and temporally linked to endosome fission. .
Tissue Specificity	
Function	similarity:Belongs to the TEX28 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night