




FABP9 rabbit pAb

Catalog No	YP-Ab-11530
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FABP9
Immunogen	Synthesized peptide derived from human FABP9 AA range: 43-93
Specificity	This antibody detects endogenous levels of FABP9 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fatty Acid Binding Protein 9;Fatty acid-binding protein 9;Testis lipid-binding protein;Testis-type fatty acid-binding protein;T-FABP;FABP9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	similarity:Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night