




CF057 rabbit pAb

Catalog No	YP-Ab-11534
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	C6orf57
Immunogen	Synthesized peptide derived from human CF057 AA range: 23-73
Specificity	This antibody detects endogenous levels of CF057 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C6orf57;SDH assembly factor 4;Succinate dehydrogenase assembly factor 4; mitochondrial;Succinate;CF057
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	similarity:Belongs to the UPF0369 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).