




APC15 rabbit pAb

Catalog No	YP-Ab-11553
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ANAPC15 C11orf51 HSPC020
Immunogen	Synthesized peptide derived from human APC15 AA range: 6-56
Specificity	This antibody detects endogenous levels of APC15 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C11orf51; HSPC020; Anaphase-promoting complex subunit 15; ANAPC15; APC15; ANAPC15 C11orf51 HSPC020
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular, nucleoplasm, anaphase-promoting complex, cytosol,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night