




FA71C rabbit pAb

Catalog No	YP-Ab-11556
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM71C
Immunogen	Synthesized peptide derived from human FA71C AA range: 38-88
Specificity	This antibody detects endogenous levels of FA71C at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA71C;FAM71C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	
Function	similarity:Belongs to the FAM71 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night