




CD019 rabbit pAb

Catalog No	YP-Ab-11566
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	C4orf19
Immunogen	Synthesized peptide derived from human CD019 AA range: 166-216
Specificity	This antibody detects endogenous levels of CD019 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD019; C4orf19
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleoplasm, cell junction,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night