




TX264 rabbit pAb

Catalog No	YP-Ab-11573
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TEX264 ZSIG11 UNQ337/PRO536
Immunogen	Synthesized peptide derived from human TX264 AA range: 112-162
Specificity	This antibody detects endogenous levels of TX264 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SIG11;TEG 264;Testis-expressed protein 264;TEX 264;UNQ337/PRO536; TX264;TEX264 ZSIG11 UNQ337/PRO536
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type III membrane protein . Cytoplasmic vesicle, autophagosome . Cytoplasm, cytosol . Nucleus . Chromosome . Is trafficked from tubular ER to growing autophagosomes via its cytosolic LIR motif (PubMed:31006537). Also found in the cytosol, nucleus and chromatin (PubMed:32152270). In response to formation of covalent DNA-protein cross-links (DPCs), localizes to the nuclear periphery, and associates with DNA replication forks (PubMed:32152270). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night