




TRI58 rabbit pAb

Catalog No	YP-Ab-11579
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRIM58
Immunogen	Synthesized peptide derived from human TRI58 AA range: 355-405
Specificity	This antibody detects endogenous levels of TRI58 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BIA2;Protein BIA2;TRIM;tripartite motif containing 58;TRI58;TRIM58
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,
Tissue Specificity	Expressed in erythroblasts.
Function	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night