




CA021 rabbit pAb

Catalog No	YP-Ab-11582
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C1orf21 PIG13
Immunogen	Synthesized peptide derived from human CA021 AA range: 24-74
Specificity	This antibody detects endogenous levels of CA021 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA021;C1orf21 PIG13
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in spleen, prostate, testis and uterus.
Function	tissue specificity:Expressed in spleen, prostate, testis and uterus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night