




F115C rabbit pAb

Catalog No	YP-Ab-11600
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM115C FAM139A
Immunogen	Synthesized peptide derived from human F115C AA range: 754-804
Specificity	This antibody detects endogenous levels of F115C at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F115C;FAM115C FAM139A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 2]: Cell membrane . Colocalizes with TRPM8 on the plasma membrane.
Tissue Specificity	Isoform 2 is expressed in the prostate and in cancerous prostate samples.
Function	similarity:Belongs to the FAM115 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night