




PHF14 rabbit pAb


Catalog No	YP-Ab-11603
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	PHF14 KIAA0783
Immunogen	Synthesized peptide derived from human PHF14 AA range: 379-429
Specificity	This antibody detects endogenous levels of PHF14 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0783;PHD finger protein 14;PHF14;PHF14 KIAA0783
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 2 PHD-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).