




CENPX rabbit pAb

Catalog No	YP-Ab-11607
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	STRA13 CENPX FAAP10 MHF2
Immunogen	Synthesized peptide derived from human CENPX AA range: 28-78
Specificity	This antibody detects endogenous levels of CENPX at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BHLHB2;BHLHE40;DEC1;SHARP 2;SHARP2;STRA13;Stra14;CENPX;STRA13 CENPX FAAP10 MHF2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome, centromere . Chromosome, centromere, kinetochore . Assembly of CENPS and CENPX and its partner subunits CENPT and CENPW at centromeres occurs through a dynamic exchange mechanism. Although exchange is continuous in the cell cycle, de novo assembly starts principally during mid-late S phase and is complete by G2. CENPX being less stably bound at the kinetochore than CENPS. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).