




EFHD2 rabbit pAb

Catalog No	YP-Ab-11611
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	EFHD2 SWS1
Immunogen	Synthesized peptide derived from human EFHD2 AA range: 13-63
Specificity	This antibody detects endogenous levels of EFHD2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EF-hand domain-containing protein D2;Swiprosin 1;Swiprosin-1;SWS1;EFHD2; EFHD2 SWS1
Calculated Molecular Weight	26kD
Cell Pathway	Membrane raft . In a mouse immature B-cell line WEHI-231. .
Tissue Specificity	Found in lymphocytes; preferentially expressed in CD8+ cells.
Function	similarity:Contains 2 EF-hand domains.,tissue specificity:Found in lymphocytes; preferentially expressed in CD8+ cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night